

# Nothing changed. Just improved.

Geistlich Bio-Oss Pen®

- ▶ Unwrap
- ▶ Moisten
- ▶ Use

For more  
information visit:  
[www.bio-oss.com](http://www.bio-oss.com)




9 out of 10 dentists are  
convinced the  
Geistlich Bio-Oss Pen®  
handles well\*

# Geistlich Bio-Oss<sup>®</sup> – The Master's Choice


## Predictable results with the No. 1 bone substitute<sup>1,2</sup>

► Additional autologous bone is not required for successful clinical results<sup>3</sup>

**NEW BONE FORMATION**



**IMPLANT SURVIVAL RATE\*\***



- An optimal scaffold for high level of osseointegration and reliable bone formation<sup>3-8</sup>
- Excellent biofunctionality<sup>9-12</sup> and aesthetic results<sup>9</sup>
- Predictable and long-term success<sup>9-16</sup>

1 Millennium Research Group, Dental Biomaterials North America, 2018 (Market research).  
 2 Millennium Research Group, Dental Biomaterials Europe, 2016 (Market research).  
 3 Jensen T et al., Clin Oral Implants Res. 2012 Mar;23(3):263-73 (Clinical study).  
 4 Lee DW et al., Int J Oral Maxillofac Implants. 2009 Jul-Aug;24(4):609-15 (Clinical study).  
 5 Orsini G et al., J Biomed Mater Res B Appl Biomater. 2005 Jul;74(1):448-57 (Clinical study).  
 6 Perelman-Karmon M et al., Int J Periodontics Restorative Dent. 2012 Aug;32(4):459-65 (Clinical study).  
 7 Scarano A et al., J Periodontol. 2004 Aug;75(8):1161-6 (Clinical study).  
 8 Traini T et al., J Periodontol. 2007 May;78(5):955-61 (in vitro after clinical study).  
 9 Buser D et al., J Periodontol. 2013 Nov;84(11):1517-27 (Clinical study).

10 Jung RE et al., Clin Oral Implants Res. 2013 Oct;24(10):1065-73 (Clinical study).  
 11 Jensen SS et al., J Periodontol. 2014 Nov;85(11):1549-56 (Clinical study).  
 12 Mordenfeld A. et al., Clin Oral Implants Res. 2010 Sep;21(9):961-70 (Clinical study).  
 13 Aghaloo TL et al., Int Journal of Maxillofac Implants 2007;22:49-70 (Clinical study).  
 14 Benic GI et al., Clin Oral Implants Res. 2009 May;20(5):507-13 (Clinical study).  
 15 Dahlin C et al., Clin Implant Dent Relat Res. 2010 Dec;12(4):263-70 (Clinical study).  
 16 Juodzbalys G et al., J Oral Rehabil. 2007 Oct;34(10):781-9 (Clinical study).

\* Geistlich Pharma AG practice test, Nov 2011  
 \*\* after 1 year

# 9 out of 10 dentists are convinced the Geistlich Bio-Oss Pen<sup>®</sup> handles well\*

The Geistlich Bio-Oss Pen<sup>®</sup> is a combination of improved handling and the proven No. 1 bone substitute material<sup>1,2</sup> Geistlich Bio-Oss<sup>®</sup>, which is known worldwide for its unequalled record of clinical success.<sup>9-16</sup>



## Very good consistency\*

89% perceived the consistency of Geistlich Bio-Oss<sup>®</sup> as “good” and “very good” after expelling the excess fluid

# The Geistlich Bio-Oss Pen® is user friendly\*


89% rated the applicator as user-friendly<sup>10</sup>

## Easy moistening\*

89% of the dentists found moistening  
the granules to be simple




# Applied in many indications, such as



Case: Dr. Patrice Imbert, France

## Lateral sinus floor elevation and larger bone defects

Geistlich Bio-Oss Pen® with large granules is recommended for larger augmentations and for sinus floor elevations.




Case: Dr. Thomas Zumstein, Switzerland


## Small defects

Geistlich Bio-Oss Pen® with small granules is recommended for smaller bony defects. Small granules ensure closer contact with the surrounding bony walls.


# Ready for use in just a few steps




1 Holding the Geistlich Bio-Oss Pen® firmly, tear off the protective seal (A) and unscrew the protection cap (B). This exposes the applicator and filter cap.




2 Do not remove the red/green filter cap from the Geistlich Bio-Oss Pen® at this step.



3 For complete wetting of Geistlich Bio-Oss®, hold the applicator with both hands and place the filter cap into the sterile and non-pyrogenic physiological saline solution or the patient blood. Completely retract and depress the plunger several times using thumb pressure until Geistlich Bio-Oss® is completely moistened.



4 Once completely wetted, gently depress the plunger to eject any excess fluid.



5 Remove the filter cap from the applicator (A). Attach the curved applicator tip by screwing it on to the end of the applicator (B). The applicator is now ready for use. Depress the plunger and apply Geistlich Bio-Oss® directly to the surgical site.

# Product range

**Geistlich**  
Biomaterials



More details about our  
distribution partners:  
[www.geistlich-biomaterials.com](http://www.geistlich-biomaterials.com)

## Manufacturer

Geistlich Pharma AG  
Business Unit Biomaterials  
Bahnhofstrasse 40  
6110 Wolhusen, Switzerland  
Phone +41 41 492 55 55  
Fax +41 41 492 56 39  
[www.geistlich-biomaterials.com](http://www.geistlich-biomaterials.com)



### Geistlich Bio-Oss Pen®

Small granules 0.25 mm – 1 mm  
Available sizes: 0.25 g ≈ 0.5 cm<sup>3</sup>, 0.5 g ≈ 1.0 cm<sup>3</sup>  
Large granules 1 mm – 2 mm  
Available sizes: 0.5 g ≈ 1.5 cm<sup>3</sup>



### Geistlich Bio-Oss®

Spongy bone substitute  
Small granules 0.25 mm – 1 mm  
Available sizes: 0.25 g ≈ 0.5 cm<sup>3</sup>, 0.5 g ≈ 1 cm<sup>3</sup>, 2 g ≈ 4 cm<sup>3</sup>



### Geistlich Bio-Oss®

Spongy bone substitute  
Large granules 1 mm – 2 mm  
Available sizes: 0.5 g ≈ 1.5 cm<sup>3</sup>, 2 g ≈ 6 cm<sup>3</sup>



### Geistlich Bio-Oss® Collagen

Spongy bone substitute/Preformed block with Collagen  
Available sizes: 50 mg (2.5 × 5.0 × 7.5 mm), 100 mg  
(5.0 × 5.0 × 7.0 mm), 250 mg (7.0 × 7.0 × 7.0 mm),  
500 mg (10.0 × 10.0 × 7.0 mm)



### Geistlich Combi-Kit Collagen

Geistlich Bio-Oss®: Collagen 100 mg  
Geistlich Bio-Gide®: 16 × 22 mm



### Geistlich Bio-Gide®

Resorbable bilayer membrane  
Small granules 0.25 mm – 1 mm  
Available sizes: 25 mm × 25 mm, 30 mm × 40 mm



### Geistlich Bio-Gide® Perio

Resorbable bilayer membrane  
with sterile templates  
Available sizes: 16 mm × 22 mm